



Caring for
your image

iMAGGO.s





IMAgo.S

Lightweight, compact and portable, the ImaGo.S is the ideal scanner for on-farm pregnancy diagnosis for bovine, swine and small ruminants.

Compatible with a standard sector probe for external use as well as a sector rectal probe for use on cattle. With a wide field of view you can diagnose quickly and with confidence.

Equipped with an intuitive touchscreen interface, the scanner is simple and easy to navigate.

Rugged and robust, with 2 interchangeable batteries, the ImaGo.S can keep up with your busy day on the farm. The comfortable carry straps offer greater comfort and freedom of movement when working with animals.

The digital conception (processor + FPGA) allows the ImaGo.S hardware and software to be upgraded in the future.

External



Rectal



Specifications

Size and Weight

Weight: 1.8kg

Screen: 6.5"

Magnesium Body

Imaging Modes

B,B+B

Zoom

Cineloop

Text and annotations

Measurements and tables

Clinical Applications

Pregnancy detection

Abdominal

Storing and Connectivity

Images and clips storage

RCA port

DC output

USB

Bluetooth

Contact us now

 www.imv-imaging.com  info@imv-imaging.com  +44 (0) 1506 460023

 facebook.com/IMVimaging  [twitter/IMVimaging](https://twitter.com/IMVimaging)