

Contact us

Unit 11
Victoria Business Park
Short Street
Southend-On-Sea
Essex, SS2 5BY



+44 (0) 1702 460 310



admin@
vescustomoptics.co.uk



Follow us



Know your endoscope

Finding the model and serial number on your scope

Finding the model and serial number on your Aohua VET endoscope



The model can be found on the control body, will start with "VET- ..."



The serial number can be found on the plug, it will start "VFV..."

Finding the model and serial number on your Aohua MVE endoscope



Model can be found on the control body, will start with "MVE- ..."




Serial number can be found on the control body where you attach the portable light source, it will start "MVE ..."

Finding the model and serial number on your scope

Finding the model and serial number on your XION Bronchoscope



The model can be found on the control body where you control the distal end of the scope.




The serial number can be found on the silver ring around the eye piece, it will have "SN" in a small box before it.

Finding the model and serial number on your Aohua processors/lightbox



The model can be found on the front of the processor at the top.



The serial number can be found at the back of the processor on the silver plate, it will be next to "Product NO: ...".

Please note

Not all model and manufacturer scopes have been incorporated in this guide. This document serves merely as a guide to facilitate the location of your model and serial number. Should you encounter any challenges in identifying either the model or serial number, we encourage you to contact the VESCO workshop at +44 (0) 1702 460 310, where our team stands ready to provide assistance. •